

# GAB2 antibody



Catalog Number: 110768

## Product name

GAB2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GAB2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human tonsillitis tissue. Positive WB detected in MCF-7 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 70-74 kDa

## Preparation

This antibody was obtained by immunization of GAB2 recombinant protein (Accession Number: NM\_080491). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110768(GAB2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:110768(GAB2 Antibody) at dilution of 1:50 (under 40x lens)