

Anti-FAM50A/XAP Rabbit antibody



Catalog Number: 168768

Product name

Anti-FAM50A/XAP Rabbit antibody

IHC: 1/100

ICC/IF: 1/100

Specificity

Human, Mouse, Rat

IP: 1/20

Antibody description

FAM50A/XAP Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FAM50A

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP


Dilutions

WB: 1/1000-1/2000

Validations



Western blot detection of FAM50A/XAP in HeLa,CHO-K1,C6,Jurkat cell lysates using FAM50A/XAP Rabbit mAb(1:500 diluted).Predicted band size:40kDa.Observed band



size:40kDa.

Western blot detection of FAM50A/XAP in Smooth Muscle,HCT116,HeLa cell lysates using FAM50A/XAP Rabbit mAb(1:500 diluted).Predicted band size:40kDa.Observed band size:40kDa.