

DDA1 antibody



Catalog Number: 109768

Product name

DDA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDA1 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, mouse liver tissue. Observed molecular weight by Western-blot: 12 kDa

Preparation

This antibody was obtained by immunization of DDA1 recombinant protein (Accession Number: NM_024050). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109768(DDA1 antibody) at dilution of 1:100



IP Result of anti-DDA1 (IP:Catalog No:109768, 4ug; Detection:Catalog No:109768 1:300) with HepG2 cells lysate 1800ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:109768(DDA1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure

DDA1 antibody



Catalog Number: 109768

Goat Anti-Rabbit IgG(H+L)