

MARK1 antibody



Catalog Number: 112478

Product name

MARK1 antibody

Immunogen

[Human MARK1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MARK1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, HeLa cells, mouse heart tissue, mouse kidney tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 85 kDa, 89 kDa, 72 kDa

Preparation

This antibody was obtained by immunization of MARK1 recombinant protein (Accession Number: NM_018650). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112478(MARK1 antibody) at dilution of 1:500



IP Result of anti-MARK1 (IP:Catalog No:112478, 4ug; Detection:Catalog No:112478 1:1000) with mouse brain tissue lysate 3440ug.