

TBP antibody



Catalog Number: 115878

Product name

TBP antibody

Immunogen

[Human TBP Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TBP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, COLO 320 cells, HeLa cells, mouse colon tissue, SKOV-3 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of TBP recombinant protein (Accession Number: NM_003194). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115878(TBP antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using TBP antibody Catalog No:115878 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).