

# XRN2 antibody

Catalog Number: 116878

---

## Product name

XRN2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

XRN2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, HeLa cells. Observed molecular weight by Western-blot: 104kd

## Preparation

This antibody was obtained by immunization of XRN2 recombinant protein (Accession Number: BC006417). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116878(XRN2 antibody) at dilution of 1:200