

# Anti-MetRS/MARS Rabbit antibody



Catalog Number: 167178

## Product name

Anti-MetRS/MARS Rabbit antibody

## Dilutions

WB: 1/1000

## Specificity

ICC/IF: 1/50

Human

FC: 1/20

## Antibody description

MetRS/MARS Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Methionyl tRNA synthetase

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

WB, ICC/IF, FC

## Validations



Western blot detection of Methionyl tRNA synthetase in K562, HeLa cell lysates using Methionyl tRNA synthetase Rabbit pAb(1:1000 diluted). Predicted band size: 101kDa. Observed band size: 101kDa.



Immunocytochemistry of Methionyl tRNA synthetase(green) in HeLa cells using Methionyl tRNA synthetase Rabbit pAb at dilution 1/50, and DAPI(blue)