

TPD52 Polyclonal Antibody



Catalog Number: 160378

Product name

TPD52 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TPD52

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human TPD52 (NP_005070.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:1000 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - TPD52 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TPD52 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.




Immunohistochemistry - TPD52 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using TPD52 antibody at dilution of 1:100 (40x lens).

TPD52 Polyclonal Antibody




Catalog Number: 160378



Immunohistochemistry - TPD52 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human tonsil using TPD52 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - TPD52 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TPD52 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - TPD52 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human placenta using TPD52 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - TPD52 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using TPD52 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TPD52 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using TPD52 antibody at dilution of 1:100 (40x lens).