

VPS26B antibody



Catalog Number: 116778

Product name

Polyclonal

VPS26B antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

VPS26B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 39kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of VPS26B recombinant protein (Accession Number: NM_052875). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116778(VPS26B antibody) at dilution of 1:1000