

# DNAJB1 antibody



Catalog Number: 110088

## Product name

DNAJB1 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

DNAJB1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IF detected in Hela cells. Positive IP detected in mouse lung tissue. Positive WB detected in mouse spleen tissue, COLO 320 cells, HEK-293 cells, HeLa cells, human liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 40kd

## Preparation

This antibody was obtained by immunization of DNAJB1 recombinant protein (Accession Number: NM\_006145). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF, IP


## Dilutions

Recommended Dilution:

## Validations




mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:110088(DNAJB1 antibody) at dilution of 1:600




Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110088(DNAJB1 antibody) at dilution of 1:100 (under 10x lens)

# DNAJB1 antibody


Catalog Number: 110088



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110088(DNAJB1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using DNAJB1 antibody Catalog No:110088 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-DNAJB1 (IP:Catalog No:110088, 4ug; Detection:Catalog No:110088 1:500) with mouse lung tissue lysate 4000ug.