

# Anti-ATM Rabbit antibody



Catalog Number: 167188

## Product name

Anti-ATM Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ATM Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human ATM

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP, ChIP

## Dilutions

WB: 1/1000


ICC/IF: 1/50

FC: 1/100

IP: 1/20

ChIP: 1/20

## Validations



Western blot detection of ATM in Hela cell lysates using ATM Rabbit pAb(1:1000 diluted).Predicted band size:351KDa.Observed band size:351KDa.