

CHRNA Polyclonal Antibody



Catalog Number: 165588

Product name

CHRNA Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CHRNA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 328-517 of human CHRNA (NP_005190.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - CHRNA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CHRNA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10001 dilution. Lysates/proteins: 26ug per lane. Blocking buffer: 4% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.