

Anti-SNF5/SMARCB1 Rabbit antibody



Catalog Number: 167288

Product name

Anti-SNF5/SMARCB1 Rabbit antibody

Specificity

Human

Antibody description

SNF5/SMARCB1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human SNF5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000-1/10000


IHC: 1/50-1/200

ICC/IF: 1/20


FC: 1/20

IP: 1/20

Validations



Immunocytochemistry of SNF5(green) in HeLa cells using SNF5 Rabbit pAb at dilution 1/50, and



DAPI(blue)

Western blot detection of SNF5 in MCF7, 3T3, HeLa cell lysates using SNF5 Rabbit pAb(1:1000 diluted). Predicted band size: 44KDa. Observed band size: 44KDa.