

Anti-TBL1 (4H2) Mouse antibody



Catalog Number: 168288

Product name

Anti-TBL1 (4H2) Mouse antibody

Dilutions

WB: 1/1000

Specificity

;ICC: 1/100

Human

Validations

Antibody description

TBL1 (4H2) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human TBL1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG2b

Applications

WB, ICC/IF



Western blot detection of TBL1X in SW480 cell lysates using TBL1X mouse mAb (1:1000 diluted). Predicted band size: 58 kDa. Observed band size: 58 kDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-TBL1X mouse mAb (dilution 1:100).