

ULK4 Polyclonal Antibody



Catalog Number: 165388

Product name

ULK4 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ULK4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ULK4 (NP_060356.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - ULK4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ULK4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure




time: 30s.

Immunohistochemistry - ULK4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human oophoroma using ULK4 antibody at dilution of 1:100 (40x lens).


ULK4 Polyclonal Antibody

Catalog Number: 165388



Immunohistochemistry - ULK4 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human kidney cancer using ULK4 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ULK4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded


human gastric using ULK4 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ULK4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using ULK4 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - ULK4 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ULK4 antibody . Blue: DAPI for nuclear staining.