

Anti-VCAM1 Rabbit antibody



Catalog Number: 167388

Product name

Anti-VCAM1 Rabbit antibody

WB: 1/2000

IHC: 1/20-1/500

Specificity

Human, Mouse, Rat

ICC/IF: 1/50

FC: 1/20

Antibody description

VCAM1 Rabbit polyclonal antibody

IP: 1/20


Validations

Preparation

Antigen: A synthetic peptide of mouse VCAM1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Immunohistochemistry of VCAM1 in paraffin-embedded Human tonsil using VCAM1 Rabbit pAb

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions



at dilution 1/100

Western blot detection of VCAM1 in 3T3 cell lysates using VCAM1 Rabbit pAb(1:1000 diluted).Predicted band size:81KDa.Observed band size:110KDa.