

IL1F8 antibody



Catalog Number: 107388

Product name

IL1F8 antibody

Specificity

Human; other species not tested.

Antibody description

IL1F8 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of IL1F8 recombinant protein (Accession Number: NM_173178). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IHC: 1:20-1:200


Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107388(IL1F8 Antibody) at dilution of 1:8000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107388(IL1F8 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107388(IL1F8 Antibody) at dilution of 1:50 (under 40x lens)