

CHAF1B Polyclonal Antibody



Catalog Number: 164788

Product name

CHAF1B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CHAF1B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 270-559 of human CHAF1B (NP_005432.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:20 - 1:100

Validations



Western blot - CHAF1B Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using CHAF1B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 30s.

Immunohistochemistry - CHAF1B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using CHAF1B antibody at dilution



of 1:200 (40x lens).


CHAF1B Polyclonal Antibody



Catalog Number: 164788

Immunohistochemistry - CHAF1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using CHAF1B antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CHAF1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat spleen using CHAF1B antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CHAF1B Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat fallopian tube using CHAF1B antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - CHAF1B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human tonsil using CHAF1B antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - CHAF1B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spleen using CHAF1B antibody at dilution of 1:100 (40x lens).