

# DNAJB7 antibody



Catalog Number: 110098

## Product name

DNAJB7 antibody

## Immunogen

[Human DNAJB7 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DNAJB7 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of DNAJB7 recombinant protein (Accession Number: NM\_145174). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110098(DNAJB7 Antibody) at dilution of 1:300