

KIF18A Polyclonal Antibody



Catalog Number: 166098

Product name

KIF18A Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to KIF18A

Preparation

Antigen: A synthetic Peptide of human KIF18A

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - KIF18A Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using KIF18A antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.