

TCEA3 antibody



Catalog Number: 115898

Product name

TCEA3 antibody

Immunogen

[Human TCEA3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TCEA3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of TCEA3 recombinant protein (Accession Number: NM_003196). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115898(TCEA3 antibody) at dilution of 1:500



IP Result of anti-TCEA3 (IP:Catalog No:115898, 4ug; Detection:Catalog No:115898 1:500) with HeLa cells lysate 520ug.