

# TFAP2A,AP-2 antibody



Catalog Number: 115998

## Product name

TFAP2A,AP-2 antibody

## Immunogen

[Human TFAP2A,AP-2 Recombinant protein \(His tag\)](#)

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

TFAP2A,AP-2 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of TFAP2A,AP-2 recombinant protein (Accession Number: NM\_001032280). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115998(TFAP2A,AP-2 antibody) at dilution of 1:300