

Anti-RAP1GAP Rabbit antibody



Catalog Number: 168698

Product name

Anti-RAP1GAP Rabbit antibody

Dilutions

1:1000-1:5000

Specificity

1:20-1:50

Human, Mouse, Rat

1:20

Antibody description

1:20

RAP1GAP Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human RAP1GAP

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG


Applications

WB, ICC/IF, FC, IP

Validations



Immunocytochemistry of RAP1GAP(green) in HeLa cells using RAP1GAP Rabbit mAb at dilution



Western blot detection of RAP1GAP in C6,3T3 cell lysates using RAP1GAP Rabbit mAb(1:1000 diluted). Predicted band size: 73kDa. Observed band size: 95kDa.