

FUT11 antibody



Catalog Number: 110798

Product name

FUT11 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FUT11 Rabbit Polyclonal antibody. Positive IHC detected in human testis. Positive WB detected in mouse lung tissue, mouse testis tissue. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of FUT11 recombinant protein (Accession Number: NM_173540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110798(FUT11 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:110798(FUT11 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:110798(FUT11 antibody) at dilution of 1:100 (under 40x lens)