

Anti-MTAP Rabbit antibody



Catalog Number: 167798

Product name

Rabbit IgG

Anti-MTAP Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

MTAP Rabbit monoclonal antibody

IHC: 1/50

Preparation

IP: 1/20

Antigen: Recombinant protein of human MTAP

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of MTAP in A549,MCF-7,U251,HL-60 cell lysates using MTAP Rabbit mAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:31kDa.