

# TFAP2D antibody



Catalog Number: 116001

## Product name

TFAP2D antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

TFAP2D Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 47-52 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_172238). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116001(TFAP2D Antibody) at dilution of 1:600

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality