

LY6G5C antibody



Catalog Number: 112401

Product name

LY6G5C antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

LY6G5C Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 16, 32, 64 kDa

Preparation

This antibody was obtained by immunization of LY6G5C recombinant protein (Accession Number: BC141637). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:112401(LY6G5C Antibody) at dilution of 1:100



IP Result of anti-LY6G5C (IP:Catalog No:112401, 4ug; Detection:Catalog No:112401 1:300) with K-562 cells lysate 1720ug.