

Anti-SLC22A17 antibody

Catalog Number: 176401

Product name

Anti-SLC22A17 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to SLC22A17

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human boct.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations




Fig1::; Western Blot analysis of extracts from K562, Jurkat cells, using BOCT Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2::; Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:100