

TRNAU1AP,SECP43 antibody



Catalog Number: 116401

Product name

TRNAU1AP,SECP43 antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

TRNAU1AP,SECP43 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis tissue. Positive WB detected in MCF7 cells, A549 cells, COLO 320 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 36 kDa,21 kDa

Preparation

This antibody was obtained by immunization of TRNAU1AP,SECP43 recombinant protein (Accession Number: NM_017846). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116401(TRNAU1AP,SECP43 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116401(TRNAU1AP,SECP43 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog

TRNAU1AP,SECP43 antibody



Catalog Number: 116401

No:116401(TRNAU1AP,SECP43 antibody) at

dilution of 1:50 (under 40x lens)