

Anti-NDUFB10 Rabbit antibody



Catalog Number: 167401

Product name

Anti-NDUFB10 Rabbit antibody

WB: 1/2000

IHC: 1/50

Specificity

Human, Mouse, Rat

ICC/IF: 1/50

FC: 1/100

IP: 1/20

Antibody description

NDUFB10 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFB10

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of NDUFB10 in K562, HeLa cell lysates using NDUFB10 Rabbit pAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 21kDa.



Immunohistochemistry of NDUFB10 in paraffin-embedded Human colon cancer tissue using NDUFB10 Rabbit pAb at dilution 1/50