

ABCD1 antibody



Catalog Number: 107801

Product name

ABCD1 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

ABCD1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human lung cancer tissue. Positive IP detected in mouse liver tissue. Positive WB detected in human brain tissue, HeLa cells, mouse liver tissue. Observed molecular weight by Western-blot: 75 kDa

Validations



Preparation

This antibody was obtained by immunization of ABCD1 recombinant protein (Accession Number: NM_000033). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107801(ABCD1 antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions


Recommended Dilution:




IP Result of anti-ABCD1 (IP:Catalog No:107801, 4ug; Detection:Catalog No:107801 1:500) with mouse liver tissue lysate 5000ug.

ABCD1 antibody

Catalog Number: 107801



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:107801(ABCD1 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:107801(ABCD1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)