

SELPLG antibody



Catalog Number: 115101

Product name

SELPLG antibody

Specificity

Human; other species not tested.

Antibody description

SELPLG Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of SELPLG recombinant protein (Accession Number: BC029782). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115101(SELPLG Antibody) at dilution of 1:600