

# EXOC1 antibody



Catalog Number: 110501

## Product name

EXOC1 antibody

## Immunogen

[Human EXOC1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EXOC1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in rat brain tissue, human brain tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 102 kDa

## Preparation

This antibody was obtained by immunization of EXOC1 recombinant protein (Accession Number: XM\_017008409). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110501(EXoc1 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110501(EXoc1 antibody) at dilution of 1:50 (under 10x lens)

# EXOC1 antibody



Catalog Number: 110501



IP Result of anti-EXoc1 (IP:Catalog No:110501, 3ug; Detection:Catalog No:110501 1:1000) with mouse brain tissue lysate 3600ug.