

Rat Syndecan-1/SDC1/CD138 (His Tag) recombinant protein



Catalog Number: 503501

General Information

Protein Construction

A DNA sequence encoding the rat SDC1 (P26260) (Met1-Lys253) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

(76+20) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 23

Molecular Mass

The recombinant rat SDC1 comprises 242 amino acids and predicts a molecular mass of 25.6 kDa. The apparent molecular mass of the recombinant protein is approximately 47 and 49 kDa in SDS-

PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat Syndecan-1 / SDC1 / CD138 Protein (His Tag)
SDS-PAGE