

SORBS1 antibody



Catalog Number: 115501

Product name

SORBS1 antibody

Immunogen

[Human SORBS1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

SORBS1 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in 3T3-L1 cells, HEK-293 cells, mouse kidney tissue, mouse liver tissue. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of SORBS1 recombinant protein (Accession Number: BC042612). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:115501(SORBS1 antibody) at dilution of 1:800



IP Result of anti-SORBS1 (IP:Catalog No:115501, 4ug; Detection:Catalog No:115501 1:800) with mouse liver tissue lysate 4800ug.