

GCDH antibody



Catalog Number: 110901

Product name

GCDH antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

GCDH Rabbit Polyclonal antibody. Positive WB detected in L02 cells, SMMC-7721 cells. Observed molecular weight by Western-blot: 44-48 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of GCDH recombinant protein (Accession Number: NM_000159). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110901(GCDH Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality