

STBD1 antibody



Catalog Number: 115701

Product name

STBD1 antibody

Immunogen

[Human STBD1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

STBD1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, A549 cells, HeLa cells, human heart tissue. Positive IP detected in A549 cells. Observed molecular weight by Western-blot: 35-43kd

Preparation

This antibody was obtained by immunization of STBD1 recombinant protein (Accession Number: NM_003943). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115701(STBD1 antibody) at dilution of 1:1000



IP Result of anti-STBD1 (IP:Catalog No:115701, 4ug; Detection:Catalog No:115701 1:1000) with A549 cells lysate 1600ug.