

Anti-Annexin V Rabbit antibody



Catalog Number: 168701

Product name

WB, ICC/IF, FC

Anti-Annexin V Rabbit antibody

Dilutions

Specificity

1:2000

Human, Mouse, Rat

1:100

Antibody description

1:50

Annexin V Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Annexin V

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of Annexin V(green) in HeLa cells using Annexin V Rabbit mAb at dilution



1/200, and DAPI(blue)

Western blot detection of Annexin V in K562,C6,3T3,HeLa cell lysates using Annexin V Rabbit mAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:36kDa.