

# Anti-PBK/SPK Rabbit antibody



Catalog Number: 168011

## Product name

Anti-PBK/SPK Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PBK/SPK Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human PBK/SPK

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/500

FC: 1/50

IP: 1/20

## Validations



Western blot detection of PBK/SPK in HeLa,A549,HL-60,U2OS cell lysates using PBK/SPK Rabbit mAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:40kDa.