

PPP3CB-specific antibody



Catalog Number: 114111

Product name

PPP3CB-specific antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

PPP3CB-specific Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive IHC detected in human heart tissue. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, A431 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 57-59 kDa

Validations



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021132). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114111(PPP3CB antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunohistochemical of paraffin-embedded human heart using Catalog No:114111(PPP3CB antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded

PPP3CB-specific antibody




Catalog Number: 114111

human heart using Catalog No:114111(PPP3CB antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-PPP3CB (IP:Catalog No:114111, 4ug; Detection:Catalog No:114111 1:500) with

mouse brain tissue lysate 4000ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed A431 cells using Catalog No:114111(PPP3CB Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)