

UBC9-Specific antibody



Catalog Number: 116511

Product name

UBC9-Specific antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBC9-Specific Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive IF detected in HepG2 cells. Positive FC detected in HepG2 cells. Positive WB detected in HEK-293 cells, human testis tissue. Observed molecular weight by Western-blot: 18kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003345).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116511(UBE2I-Specific antibody) at dilution of



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116511(UBE2I-Specific antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded


UBC9-Specific antibody

Catalog Number: 116511

human lung cancer using Catalog
No:116511(UBE2I-Specific antibody) at dilution of
1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using
Catalog No:116511(UBE2I-Specific Antibody) at
dilution of 1:25 and Rhodamine-Goat anti-Rabbit



1X10⁶ HepG2 cells were stained with 0.2ug
UBE2I-Specific antibody (Catalog No:116511, red)
and control antibody (blue). Fixed with 90%
MeOH blocked with 3% BSA (30 min). Alexa Fluor
488-conjugated AffiniPure Goat Anti-Rabbit
IgG(H+L) with dilution 1:1500.