

Anti-Tuberin Rabbit antibody



Catalog Number: 167911

Product name

Rabbit IgG

Anti-Tuberin Rabbit antibody

Applications

Specificity

WB, IP

Human, Mouse

Dilutions

Antibody description

WB: 1/1000

Tuberin Rabbit monoclonal antibody

IP: 1/20

Preparation

Antigen: Recombinant protein of human Tuberin

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Tuberin in HeLa,A549,HL-60,U2OS cell lysates using Tuberin Rabbit mAb(1:1000 diluted).Predicted band size:201kDa.Observed band size:201kDa.