

# C16orf93 antibody



Catalog Number: 108611

## Product name

C16orf93 antibody

## Immunogen

[Human C16orf93 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

C16orf93 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of C16orf93 recombinant protein (Accession Number: NM\_001014979). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108611(C16orf93 Antibody) at dilution of 1:700