

Anti-Caspase 10 Rabbit antibody



Catalog Number: 167311

Product name

Anti-Caspase 10 Rabbit antibody

WB: 1/2000

IHC: 1/20

Specificity

Human

FC: 1/50

IP: 1/20

Antibody description

Caspase 10 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Caspase-10

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC, IP

Dilutions

Validations



Immunohistochemistry of Caspase-10 in paraffin-embedded Human lung cancer tissue using Caspase-10 Rabbit pAb at dilution 1/20



Western blot detection of Caspase 10 in K562 cell lysates using Caspase 10 Rabbit pAb(1:1000 diluted). Predicted band size: 59kDa. Observed band size: 59kDa.