

ACTN3 antibody



Catalog Number: 107711

Product name

ACTN3 antibody

Immunogen

[Human ACTN3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ACTN3 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of ACTN3 recombinant protein (Accession Number: NM_001104). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107711(ACTN3 Antibody) at dilution of 1:300