

C2orf27B antibody



Catalog Number: 108711

Product name

C2orf27B antibody

Clonality

Polyclonal

Immunogen

[Human C2orf27B Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

C2orf27B Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.


Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of C2orf27B recombinant protein (Accession Number: NM_214461). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108711(C2orf27B Antibody) at dilution of 1:500