

BAP31 antibody



Catalog Number: 108421

Product name

BAP31 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000


IHC: 1:20-1:200

Antibody description

BAP31 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in HEK-293 cells, A431 cells, HeLa cells, human placenta tissue, mouse placenta tissue. Positive IF detected in HeLa cells. Positive IHC detected in human colon cancer tissue. Positive FC detected in HeLa cells. Observed molecular weight by Western-blot: 28 kDa

IF: 1:10-1:100

Validations



Preparation

This antibody was obtained by immunization of BAP31 recombinant protein (Accession Number: NM_001256447). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108421(BCAP31 antibody) at dilution of 1:2000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF, IP, FC, WB


Dilutions




Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:108421(BCAP31 antibody) at dilution of 1:100 (under 10x lens)

BAP31 antibody

Catalog Number: 108421



1X10⁶ HeLa cells were stained with 0.2ug BCAP31 antibody (Catalog No:108421, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.



Immunofluorescent analysis of HeLa cells, using BCAP31 antibody Catalog No:108421 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

IP Result of anti-BCAP31 (IP:Catalog No:108421, 3ug; Detection:Catalog No:108421 1:1000) with A431 cells lysate 3300ug.

