

SUDS3 antibody



Catalog Number: 115821

Product name

SUDS3 antibody

Immunogen

[Human SUDS3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SUDS3 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HL-60 cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of SUDS3 recombinant protein (Accession Number: BC112000). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115821(SUDS3 Antibody) at dilution of 1:600