

# Anti-CD55 Mouse antibody

Catalog Number: 169121

## Product name

Anti-CD55 Mouse antibody

## Specificity

Human

## Antibody description

CD55 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human CD55 (AA: 238-357) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

FC


## Validations



Flow cytometric analysis of K562 cells using CD55 mouse mAb (green) and negative control (red).



Western blot analysis using CD55 mAb against human CD55 recombinant protein. (Expected MW



is 38.5 kDa)

Western blot analysis using CD55 mouse mAb against Raji (1) and K562 (2) cell lysate.