

GDAP2 antibody



Catalog Number: 110921

Product name

Polyclonal

GDAP2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

GDAP2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 56-60 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of GDAP2 recombinant protein (Accession Number: NM_001135589). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:110921(GDAP2 Antibody) at dilution of 1:1000